




Description:



- ① Approx 2mm to apical foramen
- ② Approx 0.6mm to apical foramen
- ③ Apex(apical foramen)

Feature:



Multi-frequency measurement technology, the apical stenosis area is enlarged and the measurement is accurate.



The body is dexterous, the angle can be adjusted, and the measurement position can be discerned according to the sound.



Easy to measure. LED backlight display, Soft colors for a comfortable look.



High quality accessories

Continuous measurement of the root canal

Technical:

Parameters: Battery: 3.7V/320mAh

Power: $\leq 0.5W$

Display: 2.4 inch LCD screen

Host weight: 51.7g